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## System component


Figure 1. A) Artist's rendering of nano-porous biosensor principle. B) Overview of optical path and physics of signal acquisition. (Martin Latterich et al, Proteome Science, 2008. 6: 31)




Figure 2. The component of Affinity Capture- Mass Spectrometry

Silicon Kinetics provides, as an add on module, it's AC-MS interface box.

- Up to 15,000 psi
- 2, 10 port/ 2 position valves
- Adjustable collection volume
- Fully integrated into SKI Report
- USB controllable
- MS controls timing through contact closure

Inhibitors of Carbonic Anhydrase II MW 95–331 Da

A. Furosemide (MW 330.74)  
B. 4-Sulfonylbenzoic acid (MW 201.2)  
C. Benzensulfoamide (MW 157.2)  
D. Methyl sulfonamide (MW 95.1)

Figure 4. Kinetic analysis of Enzyme-small molecule

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## Results and discussion

### Immobilization and Kinetics analysis of CAII/Furosemide


Figure 5 Immobilization of CAII with EDC/sNHS on COOH Tip. Here using approximately 400 µg/mL getting 50 nm on surface. A: Sample, B: Reference

Table 1. LC/MS conditions of Affinity capture-Mass Spectrometry

Name	k <sub>off</sub> 1/s	k <sub>on</sub> 1/uM*s	K <sub>d</sub> uM
Furosemide	0.5241864	0.4523889	1.158708

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LC/MS components		Mass Spectrometry conditions	
nano-ESI LC/MS/MS	ZAPLOUS System, (AMR Inc.)	Ion Source	Dream spray(AMR Inc.)
HPLC	Paradigm MS4 (Microm BioResources, Inc.)	Mode	SRM, Full MS Scan
Autosampler	HTS-PAL(CTC Analytics)	Polarity	Negative
Loading Pump	Agilent 1100(Agilent technologies)	IS	1.6 keV
Mass Spectrometry	Thermo LTQ(ThermoFisher Scientific)		

LC conditions	
Mobile A	0.1% HCOOH/2% Acetonitrile
Mobile B	0.1% HCOOH/90% Acetonitrile
Gradient	B%: 5(0min)-95(10min)-95(12min)-5(12.5min)-5(15min)
Flow rate	5ul/min
Flow rate (loading pump)	10ul/min

Table 2. LC/MS conditions of Affinity capture-Mass Spectrometry


Figure 6. MS spectrum of Furosemide by AC-MS analysis.

A. SRM chromatogram of 0.5uM Furosemide. B. SRM chromatogram of 0 uM Furosemide. C. SRM profile of Furosemide(m/z); Precursor MS/Product MS : 329/285) D. Full MS spectrum profile, E. MS/MS spectrum of Furosemide

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- 参考文献
- 1) Martin Latterich and Jacques Corbeil, Label-free detection of biomolecular interactions in real time with a nano-porous silicon-based detection method. Proteome Science, 2008. 6: 31
  - 2) Rebecca L. Richa and David G. Myszka, Grading the commercial optical biosensor literature—Class of 2008: ‘The Mighty Binders’. J. Mol. Recognit. 2010; 23: 1–64
  - 3) Iain Campuzano, Kiaus Michelsen & Paul Schinier, On-line Nanopore Optical Interferometry Mass Spectrometry for Screening and Quantifying Small Molecule-Protein Interactions. ASMS 2013 Minneapolis, Minneapolis, MN, USA.
  - 4) Giuseppe A. Papalia , Stephanie Leavitt, et al. Comparative analysis of 10 small molecules binding to carbonicanhydrase II by diVerent investigators using Biacore technology. Analytical Biochemistry 359 (2006) 94–105